Item No	Item name
451150	EquiDress, 10 items / pk

Languages

Language – click at relevant language below		Page
English	GB • User instructions / Information	1
German	DE • Gebrauchsanleitung Warenzeichen	2
French	FR • Mode d'emploi d' EquiDress	4
Espain	SP • Modo de uso de EquiDress	6
Italy	IT • Come utilizzare Equidress	8

GB • User instructions / Information

How to use EquiDress

Before applying the dressing always clean dirt and mud from the affected area - no matter if the dressing is used as warm wet, cold wet or dry.

Warm Wet Dressing

1. Begin by heating water until it reaches a boiling point.
2. Pour the boiling water into a shallow container and allow it to cool down to 37°C.
3. Cut the dressing to the desired size and shape.
4. Submerge the dressing in the lukewarm water until it becomes saturated.
5. Once saturated, remove it from the water and gently squeeze out any excess liquid.
6. Place the dressing over the affected area, making sure the shiny(plastic) side faces away from the skin.
7. Secure the dressing in place with an appropriate bandage.

Cold Wet Dressing

1. Follow the same steps as for Warm Wet Dressing to soak the dressing in boiled water. But...
2. Allow the soaked dressing to cool completely to room temperature.
3. Apply the cooled dressing to the affected area, ensure the shiny(plastic) side faces away from the skin.
4. Optionally, you can cool the dressing in the refrigerator and store it in a bag for future use.

Dry Dressing

1. Apply the dry dressing directly to the affected area.
2. Ensure that the shiny(plastic) side of the dressing is facing away from the skin in all cases.
3. Secure the dressing with a cohesive bandage, taking care not to wrap it too tightly.

“The use of poultices in wound/injury management and other infected and non-infected sites has been commonly practiced over many centuries dating back to the days when kaolin and sheep wool poultices were the treatment of choice to "draw out" infection. A poultice is applied damp for cooling, warm for stimulation of blood flow, and dry for “drawing out pus and exudate” but must conform to the tissues it is applied over. Applied wet it also hydrates tissues such as the hoof capsule and so are used under a suitable bandage cover,

Poultices are now most commonly used as a means of “drawing out” infection from abscess sites or infected injuries and for softening the hoof capsule, They are also applied as a precautionary measure post-exercise “to prevent leg filling” and to support tendons and ligaments under a stable bandage on the lower limb regions where there is no obvious injury at all. Poultices can be used warm in areas where extra circulatory flow needs to be encouraged, Poultice materials are used to warm or cool the area under different conditions; they are a passive but valuable part of equine management both in damaged and undamaged tissues.

Historically chemicals such as boric acid and even antibacterial powders have been added to the construction of commercial poultices; this is presumably because of the supposed benefit of antiseptic effects. However, within the poultice itself these additives probably have no material benefit and, in fact, could be more harmful than helpful because this gives a false impression of its therapeutic value. A poultice is best viewed as a means of warming or cooling tissues and so acts more as a passive means of encouraging or reducing the blood flow in an area of inflammation. There is therefore no significant value in additive materials and where antibacterial properties are required other specific measures must be taken.

Professor Derek C Knottenbelt OBE, BVM&S, DVM&S, DipECEIM, DACVIM(LA), MRCVS

RCVS & European Specialist in Equine Internal Medicine.

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